



iProbe Control Software

iProbe allows three-dimensional control via two-handed operation – your mouse for X and Y movements and a jog dial for Z movements. The probe tip becomes your mouse pointer which gives you a familiar and intuitive feeling, even including dynamic movement response.



The iProbe software is integrated into the Advanced Probing Tools software suite. At the same time, it is available as a stand-alone tool.

Its streamlined display provides access to various indicators at a glance. Most importantly: which probe tip is currently active as well as the micromanipulators' speed setting. Additionally, a number of predefined macros allow fast and reliable execution of routinely required tasks.

The USB jog shuttle device yields access to the software's most important functions via five easy-to-reach buttons. The dial allows driving the currently active probe tip in the Z direction.



Further information

- Contact us at info@kleindiek.com
- Find your local agent at www.kleindiek.com